

ULTRASONIC CLEANER

MODEL: GT SONIC-D Series

FEATURES

- Normal and soft power control to resolve the cleaning blind spoit.
- Degas function for imporve the cleaning effect.
- Time Setting: 1~99 minutes adjustable.
- Heating: 0~80°C temperature adjustable
- User-friend and clear panel
- Memory function
- LED display for temperature and timer
- High-performance transducer
- Patent products

ST SURE MANUAL STANDARD CENTS The standard of the standard before press The standard of the standard of the standard before press The standard of the standard of

APPLICATIONS

GT SONIC D-Series ultrasonic cleaner are used widely in laboratory, Dental, Gem/Jewelry, Optics, Electronics, Automotive, Industry, Pharmaceutical, to provide simple, effective results for cleaning objects.

Degas function especially suitable for "beyond cleaning" applications, such as sample preparation, which includes degassing, mixing, homogenization, dissolving, cell lysing and dispersion of particles.







บริษัท โปรโทรนิกส์ อินเตอร์เทรด จำกัด

Picture	or sone	nor.	ar tour	arton.			
Model	GT SONIC-D2	GT SONIC-D3	GT SONIC-D6	GT SONIC-D9	GT SONIC-D13	GT SONIC-D20	GT SONIC-D27
Tank Size (LxWxH) mm	150×140×100	240×140×100	300×155×150	300×240×150	330×300×150	500×300×150	500×300×200
Unit Size (LxWxH) mm	190×170×220	270×170×240	330×180×310	330×270×310	330×300×150	500×300×150	500×300×200
Volume	2L	3L	6L	9L	13L	20L	27L
Ultrasonic Power	50W	100W	150W	200W	300W	400W	500W
Ultrasonic Frequency				40 kHz			
Heating Power	100W	100W	300W	300W	400W	500W	500W
Timer				1~99 Min		·	
Heating Temp.				0~80°C			

ACCESSORIES







<u>Lid</u>* Model D2 Lid is plastic